



MEDICAL, LABORATORY, AND SCIENTIFIC REFRIGERATORS AND FREEZERS

Models: EQFR Series Chest Refrigerators EQFF Series Chest Freezers VT/VLT Series Low Temperature Freezers

User Manual

BEFORE USE, PLEASE READ AND FOLLOW ALL SAFETY RULES AND OPERATING INSTRUCTIONS

Accucold Division of Felix Storch, Inc.
An ISO 9001:2015 registered company
770 Garrison Avenue
Bronx, NY 10474
www.accucold.com

TABLE OF CONTENTS

Appliance Safety	3
Important Safeguards	3-4
Installation Instructions	5-6
Locating your Appliance	5
Leveling your Appliance	5
Grounding your Appliance	6
Key Eject Lock	6
Reversing the Door Swing	6
Operating your Appliance	7
Setting the Temperature	8-10
For Dial Thermostats	8
For Digital Thermostats	8
Traceable Thermometer Alarm	9
Keypad Lock	10
Energy-Saving Tips	10
Operating Specific Models	11
Cleaning and Care	12-13
Cleaning a Frost-Free Appliance	12
Cleaning a Cycle or Manual Defrost Appliance	12
Cleaning the Base Grill	12
Changing the Light Bulb	12
Loss of Power	13
Moving your Appliance	13
Assistance or Service	13
Troubleshooting	14
Limited Warranty	16

APPLIANCE SAFETY

Your safety and the safety of others are very important.

We have provided many important safety messages in this manual and on your unit. Always read and obey all safety messages.



This is the Safety Alert Symbol. The symbol alerts you to potential hazards that can kill or injure you and others. All safety messages will follow the Safety Alert Symbol and either the words "DANGER" or "WARNING".



DANGER means that failure to heed this safety statement may result in severe personal injury or death.



WARNING means that failure to heed this safety statement may result in extensive product damage, serious personal injury, or death.

All safety messages will alert you about the potential hazard, tell you how to reduce the chance of injury, and let you know what can happen if the instructions are not followed.

IMPORTANT SAFEGUARDS



Before the equipment is used, it must be properly positioned and installed as described in this manual, so read the manual carefully. To reduce the risk of fire, electrical shock or injury when using this equipment, follow basic precautions, including the following:

⚠ DANGER **⚠**

- Plug into a grounded 3-prong outlet, do not remove grounding prong, do not use an adapter, and do not use an extension cord.
- Replace all panels before operating.
- It is recommended that a separate circuit serving only your unit be provided. Use receptacles that cannot be turned off by a switch or pull chain.
- Never clean the equipment parts with flammable fluids. These fumes can create a fire hazard or explosion. And do not store or use gasoline or other flammable vapors and liquids in the vicinity of this or any other pieces of equipment. The fumes can create a fire hazard or explosion.
- Before proceeding with cleaning and maintenance operations, make sure the power line of the unit is disconnected.
- Do not connect or disconnect the electric plug when your hands are wet.
- Unplug the unit or disconnect power before cleaning or servicing. Failure to do so can result in electrical shock or death.
- Do not attempt to repair or replace any part of your unit unless it is specifically recommended in this manual. All other servicing should be referred to a qualified technician.
- This unit is CFC- and HFC-free and contains a small quantity of Isobutane (R600a), which is environmentally friendly but flammable. It does not damage the ozone layer, nor does it increase the greenhouse effect. care must be taken during transportation and setting up of the unit that no parts of the cooling system are damaged. leaking coolant can ignite and may damage the eyes.
- In the event of any damage:
 - Avoid open flames and anything that creates a spark,
 - Disconnect from the electrical power line,

- o Air the room in which the unit isolated for several minutes, and
- Contact the Service Department for advice.
- The more coolant there is in a unit, the larger the room it should be installed in. In the event of a leakage, if the unit is in a small room, there is the danger of combustible gases building up. For every ounce of coolant, at least 325 cubic feet of room space is required. The amount of coolant in the unit is stated on the data plate inside the unit. It is hazardous for anyone other than an Authorized Service Person to carry out servicing or repairs to this piece of equipment.
- Take serious care when handling, moving, and using the unit to avoid either damaging the refrigerant tubing or increasing the risk of a leak.
- Replacing component parts and servicing shall be done by factory authorized service personnel so as to minimize the risk of possible ignition due to incorrect parts or improper service.

⚠ WARNING ⚠

FOLLOW WARNING CALLOUTS BELOW ONLY WHEN APPLICABLE TO YOUR MODE!

- Use two or more people to move and install unit. Failure to do so can result in back or other injury.
- To ensure proper ventilation for your unit, the front of the unit must be completely unobstructed.
 Choose a well-ventilated area with temperatures above 60°F (16°C) and below 90°F (32°C). [For optimal performance, install the unit where the ambient temperature is between 72° and 78°F (23°-26°C).] This unit must be installed in an area protected from the elements, such as wind, rain, water spray or drips.
- The unit should not be located next to ovens, grills or other sources of high heat.
- The unit must be installed with all electrical, water and drain connections in accordance with state and local codes. A standard electrical supply (115 V AC only, 60 Hz), properly grounded in accordance with the National Electrical Code and local codes and ordinances, is required.
- Do not kink or pinch the power supply cord of the unit.
- The size of the fuse (or circuit breaker) should be 15 amperes.
- It is important that the equipment be leveled in order to work properly. You may need to make several adjustments to level it.
- All installations must be in accordance with local plumbing code requirements.
- Make certain that the pipes are not pinched, kinked or damaged during installation.
- · Check for leaks after connection.
- Never allow children to operate, play with or crawl inside the unit.
- Do not use solvent-based cleaning agents or abrasives on the interior. These cleaners may damage or discolor the interior.
- Use this equipment only for its intended purpose as described in this *Instruction Manual*.
- Keep fingers out of the "pinch point" areas. Clearances between the door and cabinet are necessarily small. Be careful closing the door when children are in the area.

⚠ DANGER ⚠

Risk of child entrapment!

Child entrapment and suffocation are not problems of the past. Junked or abandoned appliances are still dangerous, even if they will "just sit in the garage a few days."

Before discarding your old refrigerator:

- Take off the doors
- o Leave the shelves in place so that children may not easily climb inside

SAVE THESE INSTRUCTIONS

INSTALLATION INSTRUCTIONS

Use two or more people to move and install your appliance. Failure to do so can result in back or other injury.

CAUTION: YOU MUST ALLOW THE APPLIANCE TO SIT UPRIGHT FOR 24 HOURS BEFORE TURNING ON, TO ALLOW THE COMPRESSOR LUBRICANT TO DRAIN. Failure to do so may adversely affect the performance and the lifetime of the unit.

- 1. Remove the packaging materials.
- 2. Remove any tape and glue residue from surfaces before turning on the appliance. Rub a small amount of liquid dish soap over the adhesive with your fingers. Rinse with warm water and dry.
- 3. Do not use sharp instruments, rubbing alcohol, flammable fluids, or abrasive cleaners to remove tape or glue. These products can damage the surface of your appliance.
- 4. Your appliance is heavy. When moving the appliance for cleaning or service, be sure to protect the floor. Always pull the appliance straight out when moving it. Do not wiggle or "walk" the appliance when trying to move it, because floor damage can occur.
- 5. After you remove all of the packaging materials, clean the inside of your appliance before using it. See cleaning instructions in the *Cleaning and Care* section.
- 6. Some upright models have leveling legs. Level as necessary for best operation.
- 7. If your appliance has casters, after moving into place, lock the brakes on those casters supplied with brakes. Always unlock before moving the appliance for cleaning.

Locating your Appliance

NOTE: Keep flammable materials and vapors, such as gasoline, away from appliances. Failure to do so can result in fire, explosion, or death.

- Do not install the appliance near an oven, radiator, or other heat source, nor in a location where the temperature will fall below 40°F (5°C). Do not place in direct sunlight.
- Check the strength of the floor before installing your appliance. It should hold the weight of a fully loaded appliance.
- Allow at least 1" (2.5 cm) of airspace at the top and on both sides of the appliance and 1" (2.5 cm) between the back of the appliance and the wall.
- If the appliance is going to be placed against a wall, leave extra space on the hinge side so the door can open wider than its width.

Leveling your Appliance

NOTE: If your model is equipped with front rollers, the front legs should be lowered enough to lift the rollers off of the floor to ensure the appliance does not roll forward when the door is opened. If your appliance is equipped with casters, the two front casters have brakes to prevent rolling. Keep brakes locked, and unlock for moving only.

- 1. Locate the two front leveling legs that may be taped among the packing materials. In many models these legs are already mounted into the underside of the appliance.
- 2. Screw these legs into the front holes on the bottom of the appliance. The back of the appliance rests on two fixed supports. (Built-in models also have rear level legs that can provide more precise leveling).
- 3. Adjust the front legs just enough to lift the cabinet up off the floor.
- 4. Place a leveling tool on top of the appliance, first side to side, then front to back. Turn the leg counterclockwise to raise a corner. Turn the leg clockwise to lower a corner.
- 5. Repeat Step 4 until the appliance is level.
- 6. After the appliance is fully loaded, recheck with level as in Step 4.

Grounding your Appliance

A 115 Volt (AC only) 60 Hz grounded electrical supply with a 15- or 20- amp fuse or circuit breaker is required. It is recommended that a separate circuit serving only your appliance be provided. Use outlets that cannot be turned off by a switch or pull chain. Do not use an extension cord.

Key Eject Lock

Some models have factory-installed lock to secure the contents of the unit. To lock or unlock the appliance, you must hold the key in the lock while turning. The lock is spring loaded, so the key cannot be left in the lock.


If the key is lost, contact your dealer or SUMMIT APPLIANCE for a replacement. Always give the model and serial number of your appliance when requesting a replacement key. Keys for most models can be ordered online at www.summitappliance.com/support

Reversing the Door Swing

For select units, if you find the direction of opening the door on your appliance inconvenient, you can change it. Holes on the opposite side have already been prepared in the factory.

For Select One-Door Units

- Take off the upper hinge cover from the door and remove the screws that hold the top hinge.
- Lift the hinge straight up to free the hinge pin from the socket in the top of the door.
- Lift the door up and away to free its bottom socket from the hinge pin.
- · Remove the screws to the bottom hinge.
- Remove the decorative caps from the opposite side and plug the holes from the original side with the decorative caps.
 - Move the bottom hinge to the opposite side and tighten the screws. Install any lock washers or nuts if they were used in the original installation.
 - Replace door onto the hinge pin on the bottom hinge. Place top hinge pin into door. Insert screws into the top hinge and tighten.



OPERATING YOUR APPLIANCE

In order to ensure proper temperature, you need to permit proper airflow in the appliance. Do not block any air vents. If the vents are blocked, airflow will be prevented, and temperature and moisture problems may occur. Always position contents to allow some airflow between shelves or baskets.

Allow the appliance to cool at least 6 to 8 hours before adding contents. Allow 24 hours for temperature to stabilize.

NOTE: When the appliance is first started, the compressor will run constantly until the appliance is cooled. It may take 6 to 8 hours (or longer), depending on the ambient room temperature and the number of times the appliance door is opened. After the appliance is cooled, the compressor will turn on and off as needed to keep the appliance at the proper temperature.

For optimum temperature performance, it is recommended that this unit be at least 30% full during use. If you are only using this unit to store small loads, we recommend adding some water bottles or freezer packs to ensure temperature stability.

The outside of your appliance may feel warm. This is normal. The appliance's design and main function is to remove heat from packages and air space inside the appliance. The heat is transferred to the room air, making the outside of the appliance feel warm. (In some models the condenser that gets warm is just inside the side walls, making these walls warm to the touch.)

SETTING THE TEMPERATURE

For Dial Thermostats

Your appliance may be controlled with a thermostat knob located inside the unit, towards the right-hand corner. Turn the knob clockwise from the STOP (0) position toward the higher numbers for colder temperatures.



- Higher settings on the knob (higher numbers) mean lower temperatures (colder)
 in all compartments of the appliance. The temperature inside a refrigerator may
 even drop below 32°F. Use higher thermostat settings only when required or recommended.
 When the ambient temperature is normal, we recommend a medium setting.
- Changes of the ambient temperature affect temperature in the interior of the appliance, so keep this in mind when choosing the setting of the thermostat knob.
- In STOP (0) position, the appliance does not operate (the cooling system is switched off), yet the power supply is not interrupted (the light is on when you open the refrigerator door).
- The temperatures in the refrigerator and in the freezer compartment also depend on how often you open the door.

NOTE: In all small refrigerators and freezers without interior fans, temperatures are coldest near the bottom, and occasionally near the rear. Place contents accordingly.

• For optimum temperature performance, it is recommended that this unit be at least 30% full during use. If you are only using this unit to store small loads, we recommend adding some water bottles or freezer packs to ensure temperature stability.

For Digital Thermostats

If you find that it is necessary to change the temperature setting for your digital-thermostat controlled appliance, please refer to the following instructions:



To see the SET POINT

Push and immediately release the SET key. The display will show the SET POINT value. Wait 5 seconds and the probe value will display again.

To Change the SET POINT

Push and hold the SET key for more than 3 seconds to change the set point value. The value of the set point will be displayed, and the "snowflake" LED will start blinking.

To Change the SET VALUE

Push the UP or DOWN arrow to the desired setting. The new setting will be memorized after 15 seconds, provided no other keys are pushed.

NOTE: The minimum and maximum range of your thermostat is preset and locked by us, so it cannot easily be accidentally moved between freezer and refrigerator temperatures. To reset this range lock, either call us or find information on our website under the *Troubleshooting* section.

Traceable Thermometer/Alarm

Some units are equipped with a Traceable Thermometer/ Alarm that simultaneously displays minimum, maximum, and current interior temperatures and provides a visual and audible signal when the temperature inside the unit falls outside the set temperature range.

The temperature sensor, or probe, is enclosed in a glycol-filled bottle that you insert into the location to be monitored, such as inside a refrigerator or freezer, once the unit has been programmed. The thermometer/alarm itself can be attached to the door of the unit with the magnet on the back, or it can be placed on a convenient table or shelf.

Power to the thermometer/alarm is supplied with a lithium iron disulfide battery with a 10-year shelf life. When it is inserted correctly, the unit will beep and all segments of the LCD display will be temporarily activated. (Note that alarm temperatures will have to be reset after replacing the battery.) The temperature display can be switched between Celsius (°C) and Fahrenheit (°F) by sliding the °C/°F button located on the back of the unit.

Function of Controls

- **MODE/RESET Key:** Press ONCE to reset Max/Min temperature record. When the unit is sounding an alarm, press ONCE to turn off the beeping sound.
- HIGH/LOW Key: Press ONCE to display high alarm setting. Press TWICE to display low alarm setting. Press THREE TIMES to display current temperature reading.
- ALERT ON/OFF Key: When the unit is sounding an alarm, press ONCE to turn off. However, the High-Alarm or Low-Alarm will still flash until the temperature goes back to normal or the Alarm function is turned off.

Programming the Alarms

The thermometer displays the internal temperature of the unit being monitored, as well as the maximum and minimum temperatures reached. The maximum and minimum readings reflect the highest and lowest temperatures since the battery was inserted or since the MODE/RESET key was last pressed. If the temperature rises above or falls below the set temperature range, the alarm will sound. By default, the high alarm temperature setting is 8°C (46.4°F) and the low alarm temperature setting is 2°C (35.6°F). To reset these limits, follow these instructions:

- High Alarm Temperature Setting: Press the HIGH/LOW key ONCE, then press MODE/RESET
 to enter the High Alarm temperature setting mode. Use the Up and Down arrows to adjust the
 temperature accordingly. Press MODE/RESET once again to confirm the setting.
- Low Alarm Temperature Setting: Press the HIGH/LOW key TWICE, then press MODE/RESET
 to enter the Low Alarm temperature setting mode. Use the Up and Down arrows to adjust the
 temperature accordingly. Press MODE/RESET ONCE to confirm the setting.
- By default, the low and high alarms should be audible. To turn alarms on/off:
- Audible High Alarm Setting: To make sure the high alarm is audible, press HIGH/LOW key ONCE (to activate HI-SET on the screen). Then press the ALERT ON/OFF key ONCE to activate the HIGH-ALARM beeping sound. Press again to deactivate.
- Audible Low Alarm Setting: To make sure that the low temperature alarm is audible, press the HIGH/LOW key TWICE to activate LOW-SET on the screen. Then press the ALERT ON/OFF key ONCE to activate the LOW-ALARM beeping sound. Press again to deactivate.

NOTE: For an activated high alarm and low alarm buzzer, the display should indicate HI-ALARM and LO-ALARM as shown here.

Keypad Lock (Optional on Most Models)

The factory-set default password is **0123**. If the default password is lost, press the (located inside the battery box) for 30 seconds, then wait for the warning sound, which signifies the factory-set default password has been restored.

Setting the Password

NOTE: Each password starts with a number from 0 to 9, in order. The first password starting with 0 is the default password. For example, use 1XXX, then 2XXX, etc. for 10 users.

Changing the Password

While the door is open, press the old password and press # within 10 seconds to confirm. Enter the new password, then press # to confirm.

NOTE: The first digit cannot be changed.

To Open the Door

While the door is closed, press the confirm



key, enter the password, then press



again to

Function Setting of the Automatic Lock

To set the auto-lock time, press the key (inside the battery box), enter a number (any number between 10-99), then press # to confirm. The auto-lock time can be set between 10-99 seconds after opening the door. If the time is set less than 10 seconds, the door will not lock automatically, only manually.

To Manually Lock the Door

Press the # key.

Deleting the Door-Opening Passwords

Press the key, enter default password followed by 0000, then press # and all the door-opening passwords will be deleted except for the default password.

NOTE: To access the battery box you need to remove only the two small screws on the top and bottom of that cover. Do not remove the 4 screws holding the keypad and cover to the appliance door.

Energy-Saving Tips

- The more airflow provided around an appliance, the less energy it will use.
- If your appliance requires manual defrosting, keep ice build-up to less than 1/2".
- Keep doors open a minimum of time, and shut securely.

OPERATING SPECIFIC MODELS

NOTE: Models with BI in the model number, or VT65/VLT650 models, are designed for built-in operation and need only minimal clearance on the sides and above the units. All other models require 4 inches of airflow on the sides and above the units. All models operate best when 2 inches is provided between the rear wall and the back of the unit.

EQFR Refrigerator Series / EQFF Freezer Series

These models are chest freezers or chest refrigerators, all completely frost-free with forced air cooling. No defrosting is normally needed. For best operation, allow free airflow within the interior (do not overfill). The thermostat is located on the bottom rear on the right side of the unit. Large sizes may have lid-activated lights. The lights will generate some heat, so be careful to close lids securely to deactivate.

VT, VLT Low-Temperature Freezer Models

These models are all manual defrost, upright or chest type freezers. Keep door closed for best operation. The thermostats for these models are always digital, so follow directions on Page 8.

CAUTION: VERY LOW TEMPERATURE FREEZERS CAN DAMAGE UNPROTECTED SKIN, SO ALWAYS WEAR GLOVES WHEN OPENING FREEZERS, AND NEVER DISCARD THE WARNING LABELS ON THE FRONT OF THE FREEZER.

These freezers will operate best when adequate airflow is provided on the top, sides and rear of the unit. The compressor on some models is fan-cooled, so be sure not to block the air vents on the compressor housing. The refrigerant in low temperature freezers is different from that used in household appliances, so never have unqualified technicians attempt to repair these units. These units will operate best under some load, so try to avoid operating empty or nearly empty. If needed, fill with ice-packs to at least 25% of the volume of the freezer.

Fixed Shelf

In some units, the center shelf is fixed in place and cannot be removed. This is because the glycolfilled bottle that contains the temperature sensors for your alarm and thermometer is best stored in this location. The center of your refrigerator or freezer will provide the most accurate temperature readings. For proper temperature monitoring, keep the bottles in place and do not remove this shelf.

CLEANING AND CARE

Cleaning a Frost-Free Appliance (including models starting with EQFF and EQFR)

Frost-free means that you do not have to defrost your appliance. Frost is removed by air moving and collecting any moisture and depositing it on the cooling coil. The cooling coil periodically warms up enough to melt the frost, and the water then flows down to the drain pan (located behind the base grill) where it is evaporated. A complete cleaning should be done at least once a year.

- 1. Unplug the appliance or disconnect power.
- 2. Remove all contents.
- 3. Wash the inside walls with a solution of mild detergent in warm water or 2 tbs (26 g) baking soda to 1 qt (0.95 L) warm water.
 - Do not use abrasive or harsh cleaners such as window sprays, scouring cleansers, flammable
 fluids, cleaning waxes, concentrated detergents, bleaches or cleansers containing petroleum
 products on plastic parts, interior and door liners or gaskets. Do not use paper towels,
 scouring pads, or other harsh cleaning tools. These can scratch or damage surfaces.
- 4. Rinse well and wipe dry with a clean soft cloth.
- 5. Use the appropriate cleaning method for the exterior finish.
 - Painted metal: Wash exteriors with a clean sponge or soft cloth and a mild detergent in warm water. Do not use abrasive or harsh cleaners, or cleaners designed for stainless steel. Dry thoroughly with a soft cloth. For additional protection against damage to painted metal exteriors, apply appliance wax (or auto paste wax) with a clean, soft cloth. Do not wax plastic parts.
- 6. Plug in the appliance or reconnect power.
- 7. Replace all containers.

Cleaning a Cycle or Manual Defrost Appliance

Cycle defrost models utilize a system where ice builds up on the rear wall and melts when the compressor is off. Manual defrost models continually build ice on inner walls. For CSS Models (completely stainless steel), only use cleaners especially made for stainless steel.

NOTE: Never use a metal object to remove ice from the rear wall of the appliance. If needed, use a plastic ice scraper or a hair dryer.

Clean surfaces using the same instructions as for frost-free models.

Cleaning the Base Grill (some models)

- 1. Open the appliance door and remove the base grill by pulling down and out from the top.
- 2. Replace the base grill by putting the clips in the openings in the appliance panel. Push in until the base grill snaps into place.

Changing the Light Bulb (some models)

NOTE: Not all appliance bulbs will fit your appliance. Be sure to replace the bulb with one of the same size, shape, base size and wattage. You can purchase bulbs from www.summitappliance.com/support.

- 1. Unplug the appliance, or disconnect the power.
- 2. Remove the light shield. Most models do not require tools.
- 3. Replace the light bulb with an appliance bulb of the same wattage.
- 4. Replace the light shield.
- 5. Plug in the appliance or reconnect the power.

Loss of Power

After the power is restored, unplug the appliance for 30 minutes to allow the refrigeration system to equalize. During this time, keep the door closed as much as possible.

NOTE: A full appliance stays cold longer than a partially filled one.

Moving your Appliance

When you get to the new location, refer to the *Installation Instructions* and *Temperature Selection* sections for information on installation and setting the controls.

- 1. Remove all containers.
- 2. Unplug the appliance.
- 3. Clean it thoroughly.
- 4. Tape the appliance shut and tape the electrical cord to the cabinet.

ASSISTANCE OR SERVICE

Before calling for assistance or service, please check the *Troubleshooting*. It may save you the cost of a service call. If you still need help, follow the instructions below.

When calling, please know the purchase date and the complete model and serial numbers of your appliance. This information will help us to respond better to your request.

Replacement Parts

If you need to order replacement parts, we recommend that you only use factory-specified parts. Factory-specified parts will fit right and work right because they are made with the same precision used to build all Accucold equipment.

Call customer service toll free in the USA:

1-888-4-MEDLAB Monday to Friday 9:00 a.m. - 5:00 p.m. (ET)

Order parts online at www.accucold.com/support

TROUBLESHOOTING

If you notice a problem with your appliance, try these solutions first to avoid the cost of an unnecessary service call.

Problem	Possible Solution
The appliance will not operate.	Replace the fuse or reset the circuit breaker. Recheck to see whether the unit is operating in 30 minutes. Frost-free models will regularly run a defrost cycle.
There is water in the defrost drain pan (on some models).	The water will evaporate. It is normal for water to drip into the defrost pan. When it is humid, expect that the water in the defrost pan will take longer to evaporate.
The motor seems to run too much.	Defrost and clean the appliance, or move the package so the door closes properly. The motor will run longer under warm conditions. At normal room temperatures, the motor will run about 40% to 80% of the time. Under warmer conditions, it will run even more. Adding a large number of items warms the appliance. The motor normally will run longer to cool the unit back down. Push the door firmly shut. If it will not shut all the way, see the door will not close completely later in this section. This prevents heat transfer and makes the motor work harder. Clean the condenser coils. If not, contact a qualified technician.
Temperature is too warm.	Move contents away from the air vents. This prevents cold air movement. Be aware that the appliance will warm up when the door is opened frequently. In order to keep the unit cool, try to get everything you need out at one time. Keep contents organized so they are easy to find and close the door as soon as the items are removed. Adding a large number of items warms the unit. It can take several hours for it to return to normal temperature.
The door will not close completely or is difficult to open.	Rearrange containers so that they fit more tightly and take up less space. Put the shelves and drawers back into their correct positions.

NOTES

LIMITED WARRANTY

ONE-YEAR LIMITED WARRANTY

Within the 48 contiguous United States, for one year from the date of purchase, when this appliance is operated and maintained according to instructions attached to or furnished with the product, warrantor will pay for factory-specified parts and repair labor to correct defects in materials or workmanship. Service must be provided by a designated service company. Outside the 48 states, all parts are warranted for one year from manufacturing defects. Plastic parts, shelves and cabinets are warranted to be manufactured to commercially acceptable standards, and are not covered from damage during handling or breakage.

5-YEAR COMPRESSOR WARRANTY

- 1. The compressor is covered for 5 years.
- 2. Replacement does not include labor.

ITEMS WARRANTOR WILL NOT PAY FOR:

- 1. Service calls to correct the installation of your appliance, to instruct you how to use your appliance, to replace or repair fuses or to correct wiring or plumbing.
- 2. Service calls to repair or replace appliance light bulbs or broken shelves. Consumable parts (such as filters) are excluded from warranty coverage.
- Damage resulting from accident, alteration, misuse, abuse, fire, flood, acts of God, improper installation, installation not in accordance with electrical or plumbing codes, or use of products not approved by warrantor.
- 4. Replacement parts or repair labor costs for units operated outside the United States.
- 5. Repairs to parts or systems resulting from unauthorized modifications made to the appliance.
- 6. The removal and reinstallation of your appliance if it is installed in an inaccessible location or is not installed in accordance with published installation instructions.

DISCLAIMER OF IMPLIED WARRANTIES; LIMITATION OF REMEDIES

CUSTOMER'S SOLE AND EXCLUSIVE REMEDY UNDER THIS LIMITED WARRANTY SHALL BE PRODUCT REPAIRAS PROVIDED HEREIN. IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO ONE YEAR. WARRANTOR SHALL NOT BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, OR LIMITATIONS ON THE DURATION OF IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS, SO THESE EXCLUSIONS OR LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.



WARNING: This product can expose you to chemicals including Nickel (Metallic) which is known to the State of California to cause cancer.

For more information go to www.P65Warnings.ca.gov

Note: Nickel is a component in all stainless steel and some other metal components.

Accucold Division of Felix Storch, Inc. An ISO 9001:2015 registered company 770 Garrison Avenue Bronx, NY 10474 www.accucold.com accuceld

For parts and accessory ordering, troubleshooting and helpful hints, visit: www.accucold.com/support

Revised Oct. 2020