

PURETHERM[®]
BY **accucold.**

Bathing Wipe Warming Cabinet

SKU: BCW100

Product Features

- **Tabletop or undercounter** warming cabinet
- **Double tempered glass door** for full visibility
- **User-reversible** door swing
- **Adjustable legs** for a level fit
- **Preset** to operate at 125°F
- **Overheat protection**
- **Holds up to 54** 8-piece packets of disposable bathing clothes
- **Stainless steel pocket racks** for easy storage
- **Optional** wall-mounted installation
- **Unbranded exterior** to allow for easy private labeling

Specifications

Height: 27 1/8"

Width: 16 7/8"

Depth: 15 1/4"

Capacity: 54 packets

Door: Glass

Net weight: 44 lbs.

Gross weight: 50 lbs.

Cabinet: Black

Door swing: RHD

Reversible: Yes

Approvals: ETL, ETL-C

Compliance: UL 61010-1, UL 61010-2-010, CAN/CSA-22.2

No. 61010-1 and 61010-2-010

Warranty: Limited 1-year parts & labor



**For more information,
contact us at 888.4.MEDLAB or info@accucold.com**

ACCUCOLD DIVISION FELIX STORCH, INC. ISO 9001:2015 CERTIFIED
770 Garrison Avenue, Bronx, NY 10474 USA • TEL 1.888.4.MEDLAB • FAX 844.478.8799
info@accucold.com • accucold.com